

iMPACT Installation Instructions

Follow the instructions below to install the latest version of the iMPACT software (v1.2):

- Ensure that you have admin rights to the computer
- Ensure that Microsoft .NET Framework 3.5 SP1 is installed on the computer
- If you already have a previous version of iMPACT already installed on the computer or a previous installation has failed, then first ensure that the iMPACT software is uninstalled. Select iMPACT from 'Control Panel/Add or Remove Programs' or 'Control Panel/Programs/Program and features/Uninstall a program'. If data has been collected that needs to be retained then copy the database file 'C:\iMPACT\Data\DataCaptureDB.sdf' to a safe location. Then delete the directory 'C:\iMPACT' and its contents.
- Download the file 'Impact v1-2 32-bit.exe' from DropBox to the computer. This can be found at <https://www.dropbox.com/s/mow9fkb8vwsv9zi/Impact%20v1-2%2032-bit.exe>.
- Double click on the downloaded executable file. This will extract the installation files. When the WinZip Self Extractor dialog is shown select a folder to unzip the installation files to and click on the "Unzip" button.
- Double click on the file 'setup.exe' in the folder where the installation files were unzipped to. This will run the iMPACT installation.
- By default the software installs to the directory 'C:\iMPACT' and puts a shortcut in an iMPACT folder under the Start menu.
- Clicking on the iMPACT shortcut will start the application. For testing log in using the user name "test" and password "impact" (note user name and password are all lower case). This user is assigned to a test site. Sites capturing data for the study will be sent a user name and password specifically for their site.

For technical or clinical queries relating to the data capture tool, contact impact@uhb.nhs.uk.